SonicWeld Rx[®]

Ultraschallgerät SonicWeld Rx[®]

SonicWeld Rx[®] Ultrasonic Unit

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Contents

1	Product Liability and Warranty	28
1.1	General Information	28
1.2	Scope of Delivery	28
1.3	Intended Use	29
1.3.1	Safe Usage	30
1.4	Warranty	30
1.5	User's Inspection	31
1.6	Hotline	31
2	Operation	32
2.1	Function	32
2.2	Safety Functions	32
2.3	Preparation and Installation	33
2.3.1	General Information	33
2.3.2	Sterile Operation	33
2.3.3	Steps for the Installation	34
2.3.4	Visual Check	35
2.3.5	Automatic and Manual Functional Test	35
2.4	Intraoperative Application	36
3	Reprocessing (Cleaning / Disinfection and Sterilisation)	38
3.1	Cleaning / Disinfection	38
3.1.1	Handpiece, by Machine	38
3.1.2	Handpiece, by Hand	39
3.1.3	Device, by Hand	39
3.2	Visual Check for Damage	40
3.3	Sterilisation	40
3.3.1	Handpiece	40
3.4	Preventive Maintenance	40
4	Maintenance	41
4.1	Troubleshooting	41
4.2	Replacements	42
4.2.1	Replacing the Lamp	42
4.2.2	Replacing the Sonotrode	42
5	Replacement Parts and Accessories	43
6	Technical Data	44
7	Disposal	44

1 Product Liability and Warranty

1.1 General Information

We thank you for having decided to buy a KLS Martin product. This product carries the CE mark, which means that it satisfies the essential requirements laid down in the EC Directive concerning medical devices.

We are the manufacturer of this product:

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group

Ludwigstaler Str. 132 · D-78532 Tuttlingen · Germany

Postfach 60 · D-78501 Tuttlingen · Germany

Tel. +49 7461 706-0 · Fax +49 7461 706-193

info@klsmartin.com · www.klsmartin.com

1.2 Scope of Delivery

Description	Item no.
SonicWeld Rx® System, complete	52-500-00-04
Case with foam padding	52-502-02-04
Ultrasonic unit	52-500-01-04
Standard sonotrode, straight	52-501-01-04
Foot switch round (activation in 360° area)	52-500-02-04
Handpiece	52-500-03-04
Flat spanner	52-502-01-04
Manual	90-267-58-32

Optional equipment:

Description	Item no.
Foot switch rectangular (reduced activation resistance)	52-500-04-04
Hand activation complete	52-504-00-04

1.3 Intended Use

The SonicWeld Rx[®] is a system to fix resorbable SonicPins Rx[®] in bones. Further information on SonicPins Rx[®] is available in the manual for the osteosynthesis material Resorb x[®].

Possible areas of application include:

Non-load bearing parts of the cranium, including bone augmentation.

A load bearing part of the cranium is, for example, the mandible.

The ultrasonic unit is a medical device class IIa in compliance with EC Directive 93/42/EEC and the applicable, national legal requirements.

- According to these requirements the device may only be used for the described application and only by qualified operators. Furthermore, the following must be observed:
 - The applicable industrial safety regulations
 - The applicable accident prevention measures
 - These instructions for use
- According to these requirements operators of the device must
 - only use materials and equipment free of defects
 - observe the intended use
 - protect themselves, the patient and others from any dangers

Sonotrodes

- The standard sonotrodes are used to apply the SonicPins Rx[®]. They are available in straight and angled versions (item 19 and 20, foldout of the front page). Angled sonotrodes facilitate work in areas that are difficult to access, for example in the posterior molar regions.
- The smoothing sonotrodes are used for smoothing Resorb x[®] implants. They are available in straight and angled versions (item 21 and 22, foldout of the front page).

1.3.1 Safe Usage

- Before use, check device for functionality and proper condition.
- In order to avoid damage due to improper installation or use and not to compromise the warranty or guarantee, the following points must be observed:
 - Use device only as described in these instructions
 - Observe information for safety and maintenance
 - Observe user instructions in accordance to all applicable standards
 - Device and equipment must only be operated by persons that have the required training, expertise and experience.
 - Keep user instructions in a place where OR personnel can access them at any time.

1.4 Warranty

Our Standard Terms and Conditions of Sale effective at the time shall apply. Agreements diverging from these Standard Terms and Conditions do not restrict the legal rights of the buyer.

Any warranty exceeding the above provisions shall require a contractual form and shall exclude component-related vandalism, software updates and consumables.

Important Notices

The product may only be repaired by Gebrüder Martin or a qualified person or firm expressly authorized by Gebrüder Martin to perform such work.

If the repair is carried out by a person or firm specially authorized by Gebrüder Martin, the operator of the product is required to obtain from the repairer a certificate with details about the nature and scope of the repair work done. This certificate must be dated and signed and include the firm's details.

In all cases where a party other than the product manufacturer performed the work, repaired products must be additionally marked with the repairer's ID label.

Improper interventions or alterations performed by third parties during the period of limitation shall void any and all warranty claims. Unauthorized actions performed on the product shall invalidate any liability claims vis-à-vis Gebrüder Martin.

If medical equipment is modified, warranty claims and any existing authorisations will lapse.

1.5 User's Inspection

Immediately upon receipt, the goods must be checked for completeness and potential damage in transit. Notice of any such damage must be given immediately.

1.6 Hotline

- Should you have any questions on how to handle the unit/product or question on its clinical application, please do not hesitate to contact the Product Management
Tel: +49 7461 706-216
Fax: +49 7461 706-350
- Should you have any technical questions, please do not hesitate to turn to our Martin Service Center:
Tel: +49 7461 706-343
Fax: +49 7461 706-203
E-Mail: msc@klsmartin.com

Our service representatives are available Monday through Friday from 8 a.m. to 5 p.m. (CET).

Should you have any questions concerning maintenance contracts or training courses, please contact our Technical Service Manager at +49 74 61 706-332, or send an e-mail to: msc@klsmartin.com.

NOTICE

To answer your technical questions as efficiently as possible, our service technicians require the serial number of the product. Therefore, please have this number at hand when contacting our hotline. It is part of the information provided on the rating plate; see item 28, foldout of the back page.

2 Operation

2.1 Function

The ultrasonic unit of the SonicWeld Rx® system converts electrical energy into mechanical oscillations (ultrasound) which cause heating on the border areas between bone and SonicPin Rx® and by doing so allow the SonicPins Rx® to be inserted.

2.2 Safety Functions

- The system can only transfer energy if the resonance frequency of the sonotrode (item 19 – 22, foldout front page) is not modified due to too high mechanical load. In case of too high mechanical load (excessive pressure on the SonicPin Rx® with the handpiece) the system will be turned off automatically.
- The system is protected against overload by built-in fuses.
- Self-test: After turning on the device the system performs an automatic function check of the complete system. After the successful check the device switches to operating mode. In this check the following points are tested:
 - Whether the connected components are complete
 - Whether the handpiece is functional
 - The basic insulation
 - Whether the internal operating software has loaded properly
 - Whether the fuses are functioning properly

▲WARNING**Life hazard for the patient due to improper power connection!**

- The device must not be combined with other devices.
 - The device must not be connected to a mains that also supplies life-supporting systems.
-

2.3 Preparation and Installation

2.3.1 General Information

Gebrüder Martin only guarantees the proper functioning of the device when the following rules are observed.

When setting-up and operating the unit, the following regulations must be complied with:

- National requirements for installation and operation
- National requirements for fire and explosion prevention

The device may only be connected to a grounded outlet with alternating current of 100 – 240 V. Adjustment of the device to different mains voltages is not required.

2.3.2 Sterile Operation

Foot switch and ultrasonic unit must be set up and operated in the unsterile area of the operation area. Handpiece, hand activation (optional) and connection cable are in the sterile area of the operating area and must therefore only be used in a sterilised state.

- Ensure that the handpiece is dry.
- Ensure that the handpiece and sonotrode have been sterilised properly, see chapter 3 "Reprocessing (Cleaning / Disinfection and Sterilisation)", page 38.

▲WARNING

Risk of infection due to unsterile use!

Improper sterilisation and unsterile use of the SonicWeld Rx® system can lead to serious health risks for the patient.

Sterilisation must be done using a validated steam sterilisation method, e. g. according to EN 285:2009 in a steriliser and validated according to ISO 17665-1:2006.

2.3.3 Steps for the Installation

CAUTION**Danger of injuries due to improper ultrasonic therapy!**

Interferences with other devices for ultrasonic diagnostics and ultrasonic therapy, respectively, are currently not known, can, however, not be excluded entirely. If necessary, appropriate measures must be taken.

- Install the ultrasonic device on even, dry and solid surfaces only.
- Connect foot switch / hand activation to the connection socket for the foot switch / hand activation (item 9, foldout front page).
- Take off protection cap from the connection cable and connect the cable to the handpiece (item 1, foldout front page).
- Use a flat spanner to check if sonotrode sits securely on the handpiece (max. tightening torque 0.03 Nm).
- Connect power supply cable to the device (item 10, foldout front page) and then to the grounded mains socket.

Hand activation (optional)

There is an optional hand activation set. Please consider that the foot switch is no longer required if you install the hand activation.

The hand activation set comes with a separate user manual. Please observe the safety and maintenance information before using the set.

NOTICE

- Ensure that no objects or other devices are placed on top of the ultrasonic unit. The ventilation openings (item 8, foldout front page) must not be covered.
 - Observe all established rules for the connection of sterilised device components to unsterilised devices!
-

2.3.4 Visual Check

- Check visually if all components used are free of defects. Damaged components must not be used and must be replaced.

NOTICE

The tip of the sonotrode must not show any signs of mechanical damage. It must not be bent or rammed (typical damage after the handpiece has fallen down).

2.3.5 Automatic and Manual Functional Test

After turning on the device by pressing the mains switch (item 12, foldout front page) the device performs an automatic functional check. After the successful check the device switches to operating mode.

In case of an assumed malfunction you can perform another, manual functional check by pressing the button **Test** (item 7, foldout front page) with the device switched on.

⚠ CAUTION

Danger of injury by burns!

Contact with the sonotrode during the check can cause burns to tissue.

Avoid contact with skin, eyes, etc!

- Press the button **Test** (item 7, foldout front page) and wait until the test sequence has finished. During this process the tip of the sonotrode must not touch any objects.
 - Display **Ready**: There are no malfunctions.
 - Display **Service**: Perform steps described in chapter 4.1 "Troubleshooting", page 41.

2.4 Intraoperative Application

CAUTION

Danger of injury by burns!

When foot switch / hand activation is activated, contact of tissue with the sonotrode can cause burns.

Avoid contact with skin, eyes, etc!

- Turn on the device on the mains switch (item 12, foldout front page).
- Wait until functional check has finished.
- After determining the drilling hole depth and diameter according to indication choose the proper SonicPin Rx[®] and drill so that the drill can be inserted up to the STOP.
Drill geometry and the design of the SonicPin Rx[®] are exactly matched. The tip of the drill is designed in a way that it can drill through the Resorb x[®] plate and the Resorb x[®] mesh. The geometries of the drills are matched to the diameter of these implants.
It is recommended to use a drill that is at least 1 mm > the intended pin length.

After turning on the device it is set to the default power setting 2.

This power setting guarantees that the SonicPins Rx[®] can be inserted ideally. If required, power settings 1, 3 and 4 can be chosen depending on bone structure, ping geometry, and the operator's working method (see chapter 6 "Technical Data", page 44).

Individual adjustments are possible in manual operation mode:

- Press button  (Manual operation, item 4, foldout front page).
- Use the buttons  /  (item 5 / 6, foldout front page) to adjust the desired power setting of the sonotrode.
The chosen power setting of the sonotrode is displayed as a bargraph next to "I".
The power setting adjusted in that way is utilized when the foot switch / hand activation is activated.
- Press the button  again or one of the buttons 1 to 4 to end the manual operation mode.
- Hold the handpiece like a pencil.


- To take out the SonicPins Rx[®] from the clip magazine: Put sonotrode tips on the SonicPin Rx[®] and take out SonicPin Rx[®]. As the SonicPins Rx[®] hold magnetically on the tip of the sonotrode they can easily be taken out of the clip magazine.

- Put the SonicPin Rx[®] fixed on the sonotrode tip in the drill hole of the plate or the mesh, respectively.
- Apply some pressure to the SonicPin Rx[®] in an axial direction to the drill hole.

- Press foot switch / hand activation only so long until the SonicPin Rx[®] has been inserted into the drill hole in the desired depth by applying light pressure. After inserting the SonicPins Rx[®] the tip of the sonotrode should not be moved for at least 3 seconds.
- Loosen sonotrode tip after the SonicPin Rx[®] has solidified by turning it in and axial way and then remove it.

- Press button  again.

3 Reprocessing (Cleaning / Disinfection and Sterilisation)

The warranty for the cleaning and sterilisation of the handpiece is limited to 250 reprocessing and sterilisation cycles.

The recommendations of the Working Group Instrument Reprocessing (AKI) provide further information on how to properly reprocess reusable medical devices:

“Instrument reprocessing - how to do it right”, see website: www.a-k-i.org.

Further information on the reprocessing of medical devices can be found on: www.klsmartin.com/Informationen_zur_Au.239.0.html.

3.1 Cleaning / Disinfection

3.1.1 Handpiece, by Machine

The handpiece of the SonicWeld Rx® can be reprocessed by machine / disinfected thermally. It can be reprocessed with programme sequences that are authorised for surgical instruments and implants. For the cleaning process observe the instructions by the manufacturers of the cleaning machines and by the manufacturers of the cleaning agents. The process including loading must guarantee sufficient removal of residuals.

Only mildly alkaline cleaning agents are authorised for use.

- Disconnect handpiece from the device and secure with protection cap.
- Wind up connection cable together with handpiece in a circle. The minimum diameter must be 20 cm. We recommend using our storage rack.
- Place the handpiece in the instrument washer in a way that the opening with the sonotrode on top of the handpiece points downwards.

This is to ensure that no cleaning fluids can collect in the handpiece.

3.1.2 Handpiece, by Hand

Clean the handpiece with a clean, lint-free cloth that has been dampened with a standard disinfection agent on basis of ethanol (50/50) or isopropanol (70/30) that is authorised for the cleaning of medical instruments.

NOTICE

Risk of damage due to improper operation!

The handpiece must not be immersed in a disinfection bath or an ultrasonic bath.

Do not use any acetone-containing disinfection agents for cleaning and disinfecting the handpiece.

-
- Disconnect handpiece from the device and secure with protection cap.

3.1.3 Device, by Hand

Cleaning of the SonicWeld Rx[®] incl. foot switch / optional hand activation is done with a clean, lint-free cloth that has been dampened with a standard disinfection agent on basis of ethanol or methanol, and that is authorised for the cleaning of medical instruments.

⚠ WARNING

Life hazard due to electric shock!

It is essential to prevent that liquids penetrate into the device!

If liquids have penetrated into the device, disconnect the device from the mains immediately.

3.2 Visual Check for Damage

Perform a visual check of all components before the system is used again and before the handpiece is sterilised.

▲WARNING

Risk of infections and injuries due to damaged components!

Damaged components must not be used!

3.3 Sterilisation

Just like the surgical instruments the sonotrode must be reprocessed. Before use, a new sonotrode in sterile packing must be used.

3.3.1 Handpiece

▲WARNING

Risk of infection due to unsterile use!

Improper sterilisation and unsterile use of the SonicWeld Rx® system can lead to serious health risks for the patient.

The handpiece can be sterilised with steam.

Sterilisation is to be done using a validated steam sterilisation method, e. g. according to EN 285:2009 in a steriliser ANSI / AAMI ST46-1993 and validated according to ISO 17665-1:2006. When using the fractionated vacuum process, sterilisation must be carried out using the 134 °C / 2 bar programme with a minimum process time of 5 minutes.

We recommend using our storage rack.

NOTICE

The operator / user is responsible for cleaning, disinfecting and sterilising the device. It is essential that national regulations and restrictions are observed.

3.4 Preventive Maintenance

The ultrasonic device including handpiece and foot switch can be operated without maintenance.

Safety-related checks are limited to visual checks of the intactness of the connection cable and the handpiece, the foot switch and the mains connection.

4 Maintenance

4.1 Troubleshooting

Malfunction	Cause	Detection	Correction
Display Service after pressing the mains switch	Loose sonotrode	Display Service	Tighten sonotrode using a flat spanner (item 15, foldout front page)
	Components not connected		Connect components
	Components defective		Replace defective components
"Screeching" sounds during the self-test or during the operation of the SonicWeld Rx®.	Residual moisture or condensate on the contact points.	Clearly audible	The sound does not affect the proper functioning of the SonicWeld Rx®. It will usually stop by itself. Maybe longer drying periods after the sterilisation should be considered.
Tip of sonotrode cannot be removed from the SonicPin Rx®	Bent or rammed sonotrode tip	Visual check	Replace the sonotrode tip
No magnetic hold of the SonicPins Rx® on the sonotrode tip			
Sonotrode tip "sticks" on the SonicPin Rx®	Too long contact time between sonotrode and pin in regard to set energy level	Pin is generally deformed, sonotrode cannot be removed residue-free	Observe cooling time. Loosen adhesion by turning sonotrode in an axial way (see chapter 2.4, page 36)
	Bent sonotrode tip	Visual check	Replace the sonotrode tip
No light on the handpiece	Lamp damaged, condensate on contact points	Visual check	Consider replacing the lamp. If the lamp on the sonotrode is damaged, the ultrasonic function will still work properly
No display	No mains connection	Visual check	Connect to mains and turn on device
	Defective fuse	Electrical test	Replace fuse
	Device is not on	Mains switch is set to "0"	Set mains switch to "I"
Display ERROR during operation	Device is overloaded (jamming / jerking of the sonotrode tip, applied pressure too high)	Energy transfer is interrupted	Do not use the handpiece for a few moments, and then press foot switch / hand activation again

4.2 Replacements

4.2.1 Replacing the Lamp

- Unscrew the cap nut on the handpiece (item 14, foldout front page) by turning it counter-clockwise.
- Pull handpiece (item 14, foldout front page) from connection cable.
- Insert new lamp in the lamp socket and put handpiece (item 14, foldout front page) back on.
- Unscrew the cap nut on the handpiece (item 14, foldout front page) by turning it clockwise.

A damaged lamp does not have a negative effect on the general functioning of the SonicWeld Rx®.

4.2.2 Replacing the Sonotrode

- Unscrew sonotrode (item 19 – 22, foldout front page) from the handpiece (item 14, foldout front page) by turning it counter-clockwise using a flat spanner, and tighten by turning it clockwise with a maximum of 0.03 Nm.

5 Replacement Parts and Accessories

Description	Item no.
Mains cable	
- Mains cable Germany, France, Austria, Scandinavia	08-024-00-39
- Mains cable Great Britain, Ireland, UAE	08-024-00-40
- Mains cable Italy	08-024-00-30
- Mains cable Japan	08-024-00-31
- Mains cable Canada, Korea, Taiwan, USA	08-024-00-29
- Mains cable Switzerland	08-024-00-28
Handpiece	52-500-03-04
Standard sonotrode, straight	52-501-01-04
Standard sonotrode, angled	52-501-02-04
Smoothing sonotrode, straight	52-501-03-04
Smoothing sonotrode, angled	52-501-04-04
Flat spanner	52-502-01-04
Lamp	52-503-03-04
Foot switch round (surgery)	52-500-02-04
Foot switch rectangular (dental)	52-500-04-04
Manual	90-267-58-32
Storage rack for the sterilisation and transport (storage rack)	55-969-38-04
Fuse (T1.25H)	52-503-04-04
Hand activation (optional)	52-504-00-04

6 Technical Data

Description	Data
Fuse value	T1.25H (T = slow-blow, H = high breaking capacity)
Frequency mains connection	50 – 60 Hz
Voltage mains connection	100 – 240 V alternating voltage with protective ground contact
Consumption current	0.6 A
Power setting 1	13 W
Power setting 2	17 W
Power setting 3	20 W
Power setting 4	33 W
Operating frequency	29.5 ± 0.5 kHz
Dimensions: Length x Width x Height	26 x 26 x 5 cm
Length of handpiece connection cable	2.95 m
Length of connection cable / foot switch	2.95 m
Number of reprocessing cycles for the handpiece	250
Operating mode	Continuous operation with intermittent loading: 2 s operation 15 s break 5 repetitions 4 min break
Environmental conditions for transport and storage	
- ambient temperature	-20 °C up to +45 °C
- Relative humidity	10 – 100 %
- Air pressure	500 – 1060 hPa
Environmental conditions for operation	
- ambient temperature	+19 °C up to +35 °C
- Relative humidity	30 – 75 %
- Air pressure	700 – 1060 hPa

7 Disposal

When the device was designed, it was tried not to use composite materials wherever possible. Thanks to this design concept the device can be easily recycled at the end of its life time. We also offer you to take back the device and to dispose of it properly

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Stuckenbrock Medizintechnik GmbH
Martin Nederland/Marned B.V.
1270 AG Huizen · The Netherlands
Rudolf Buck GmbH
KLS Martin UK Ltd.