Ultraschallgerät SonicWeld Rx®

Ultrasonic Unit SonicWeld Rx®

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Datum:

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Contents

Terms & Acronyms	45
1 Product Liability and Warranty	45
1.1 General Information	45
1.2 Scope of Delivery	45
1.3 Intended Use.....	46
1.4 Warranty	47
1.5 User's Inspection	47
1.6 Hotline.....	48
2 Notices Concerning this Document.....	48
2.1 Symbols Used in this Document	49
3 Users	49
4 Operation	50
4.1 Function	50
4.2 Safety Functions	50
4.3 Preparation and Installation.....	51
4.3.1 General Information	51
4.3.2 Setting up SonicWeld Rx®	51
4.3.3 Changing the Display Language	53
4.3.4 Visual Check.....	54
4.3.5 Automatic Functional Test	54
4.4 Intraoperative Use	55
4.4.1 Sterile Working	55
4.4.2 Preparatory Measures	56
4.4.3 Application	57
4.5 Putting the Unit out of Operation.....	59
5 Processing, Cleaning, Disinfection and Sterilization	60
5.1 Cleaning/Disinfecting	60
5.1.1 Handpiece, by Machine	60
5.1.2 Handpiece, by Hand.....	61
5.1.3 Device, by Hand.....	61
5.2 Visual Check for Damage	62
5.3 Sterilization	62
5.4 Maintenance.....	63

6	Maintenance	64
6.1	Troubleshooting	64
6.1.1	Displays	64
6.1.2	Handpiece Statuses	65
6.1.3	Potential Malfunctions	66
6.2	Replacements	67
6.2.1	Replacing the Sonotrode	67
6.2.2	Replacing the Mains Fuse	67
6.2.3	Replacing the Activation LED	67
6.2.4	Replacing the Mains Cable	67
6.2.5	Replacing the Handpiece with Connecting Cable	68
7	Accessories & Spare Parts	69
8	Technical Data	70
9	Guidelines and Manufacturer's Declaration Regarding Electromagnetic Compatibility (EMC)	71
10	Disposal of the Device.....	74
11	Safety Checks (SC)	75

Terms & Acronyms

Designation	Description
EMC	Electromagnetic compatibility
LED	Light-emitting diode
REF	Reference number for ordering Gebrüder Martin products (item number)
SC	Safety Check
SonicPin	Designates a resorbable SonicPin Rx® or SonicPin xG

1 Product Liability and Warranty

1.1 General Information

We thank you for having decided to buy one of our products. This product carries the CE mark, which means that it satisfies the essential requirements laid down in the EC Directive concerning medical devices.

We are the manufacturer of this product:

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · D-78532 Tuttlingen · Germany

Postfach 60 · D-78501 Tuttlingen · Germany

Tel. +49 7461 706-0 · Fax +49 7461 706-193

info@klsmartin.com · www.klsmartin.com

1.2 Scope of Delivery

Designation	REF
SonicWeld Rx® system, comprising:	52-500-20-04
Case with foam plastic insert for SonicWeld Rx®	
SonicWeld Rx® ultrasonic unit	
1 handpiece with connecting cable (2.95 m)	
Standard sonotrode, straight	
Open-end wrench for sonotrodes	
Instructions for Use, SonicWeld Rx®	

1.3 Intended Use

The SonicWeld Rx® system is used for fixing resorbable SonicPins in the bone and for smoothing resorbable Resorb x® or Resorb xG foils and membranes.

For additional information on SonicPins as well as Resorb x® and Resorb xG implants, please refer to the "Osteosynthesis Implants Resorb x® and SonicPin Rx" Instructions for Use.

According to the EC Directive 93/42/EEC and applicable national legal provisions, the following applies:

- The system may only be used by competent medical staff for the application described, paying due regard to
 - current safety-at-work regulations
 - current accident prevention regulations
 - these Instructions for Use
- According to the above provisions, it is the user's responsibility to
 - ensure that all working equipment used is in perfect condition
 - protect himself/herself, patients and third-parties from hazards
- Gebrüder Martin shall not be responsible for complications caused by wrong indications, wrong implant selection, incorrect combination of system components and surgical technique, limitations of the treatment used, or lack of hygiene.
- The product must be used according to indication.
- The operational safety of the device must be verified at regular intervals, see section 6 "Maintenance", page 64 and section 11 "Safety Checks (SC)", page 75.
If the device is not functionally reliable and/or safe to operate, it must be marked as such and withdrawn from service. A technical check is mandatory in any such case.
- In order to avoid damage due to improper installation or use and not to compromise the warranty or guarantee, this Gebrüder Martin system must not be used in conjunction with third-party systems.

Risk of injury due to unauthorized modification of the product!

Altering the product in any way can pose a life hazard due to electric shock! Therefore, unauthorized modifications are strictly forbidden under any circumstances.

Sonotrodes

- The standard sonotrodes are used for the application of SonicPins. They are available in 2 different designs, straight and angled (11, 13, front fold-out page). The angled sonotrodes facilitate working in regions that are hard to access, e.g. the posterior teeth region.
- The smoothing sonotrodes are used for smoothing out Resorb x® and Resorb xG implants. They too are available in a straight version and an angled one (14, 15, front fold-out page).
- The angled standard sonotrode (13, front fold-out page) has an axial working end for inserting the SonicPins, and an angled working end for smoothing out Resorb x® and Resorb xG implants.

Sonotrode materials coming into body contact with the patient or user during application all have good biological compatibility in accordance with ISO 10993-1.

1.4 Warranty

Our Standard Terms and Conditions of Sale effective at the time shall apply. Agreements diverging from these Standard Terms and Conditions do not restrict the legal rights of the buyer.

Any warranty exceeding the above provisions shall require a contractual form and shall exclude component-related vandalism, software updates and consumables.

Important Notices

The product may only be repaired by Gebrüder Martin or a qualified person or firm expressly authorized by Gebrüder Martin to perform such work.

If the repair is carried out by a person or firm specially authorized by Gebrüder Martin, the operator of the product is required to obtain from the repairer a certificate with details about the nature and scope of the repair work done. This certificate must show the date of the repair and the details of the person or firm carrying out the work and must be signed.

In all cases where a party other than the product manufacturer performed the work, repaired products must be additionally marked with the repairer's ID label.

Improper interventions or alterations performed by third parties during the period of limitation shall void any and all warranty claims. Performing unauthorized actions on the product is strictly prohibited at any time. Non-compliance will void any liability claim against Gebrüder Martin.

1.5 User's Inspection

- Immediately upon receipt, the goods must be checked for completeness and potential damage in transit.
- Notice of any such damage must be given immediately.

1.6 Hotline

- Should you have any questions on how to handle the product or use it for clinical applications, please do not hesitate to contact the Product Management:

Tel: +49 7461 706-216

Fax: +49 7461 706-350

NOTICE

To answer your technical questions as efficiently as possible, our service technicians require the serial number of the product. Therefore, please have this number at hand when contacting our hotline. It is part of the information provided on the rating plate (16, rear fold-out page).

2 Notices Concerning this Document

Non-observance of this document can lead to serious or even lethal patient or user injury!

Improper handling and care as well as non-intended use can lead to premature wear and/or pose a risk to patients and users!

It is the operator's responsibility to ensure that all personnel handling the product have understood and do observe the notes and instructions in this document.

- Every user is required to read this document completely and follow it carefully.
- In particular, be sure to heed all cautions, warnings and danger notices.
- Keep this document accessible to users at all times.
- To ensure the functional reliability of the SonicWeld Rx®, the "Osteosynthesis Implants Resorb x® and SonicPin Rx" Instructions for Use must be observed as well.
- It is the user's responsibility to inform the patient on all the application, cautionary, warning and danger notices contained in this document and to ensure that he/she has understood them.

2.1 Symbols Used in this Document

Throughout this document, important information (such as general or safety-related notices) is marked with the following symbols and signal words:

▲WARNING

Life hazard or danger of serious injury!

Indicates a situation which, if not avoided, **could** result in death or serious injury!

▲CAUTION

Risk of injury!

Indicates a situation which, if not avoided, could result in minor or moderate injury!

NOTICE

Risk of material damage!

Indicates a situation which, if not avoided, could lead to material damage (loss of time, data loss, device/machine failure, etc.)!

3 Users

The SonicWeld Rx® system may only be used by competent medical staff who, due to their specialized training, experience and knowledge of the relevant regulations, have the necessary qualification for using the system correctly.

Users must have adequate knowledge of the rules of medical practice, as well as the state of the art and the contents of pertinent scientific articles published by medical authors.

The user must be thoroughly familiar with the surgical technique as well as with the set of implants and the instruments and accessories used.

The user bears the responsibility for

- the proper performance of the surgical procedure
- selecting the right implant components and implanting them properly
- avoiding or minimizing the general risks posed by surgical interventions

4 Operation

4.1 Function

The ultrasonic unit of the SonicWeld Rx® system converts electric energy into mechanical vibrations (ultrasound). When using a standard sonotrode, the ultrasonic energy heats up the resorbable material at the interfaces between the bone and the SonicPins, thus causing the SonicPins to glide into the predrilled hole. When using a smoothing sonotrode, the ultrasonic energy heats up the Resorb x® or Resorb xG mesh, plate, foil or membrane implants, thus allowing to smooth them out.

4.2 Safety Functions

Energy transmission is possible only if the resonance frequency of the sonotrode (11, front fold-out page) is not altered by excessive mechanical load (e.g. when applying too much pressure to the SonicPin with the handpiece). In such case of overload, the system is switched off automatically.

Built-in mains fuses protect the unit from overload.

The device automatically performs a self-test during start-up, see section 4.3.5 "Automatic Functional Test", page 54.

During the self-test, the system checks that:

- all components required are connected
- the handpiece functions properly
- the basic insulation is intact
- the internal operating software has been loaded correctly
- the mains fuses are intact

Life hazard for the patient due to improper power connection!

- The device must not be combined with other devices.
 - The device must not be connected to a mains that powers life-supporting systems.
-

4.3 Preparation and Installation

4.3.1 General Information

Gebrüder Martin only guarantees the proper functioning of the device only if the following rules are observed:

- National requirements for installation and operation
- National requirements for fire and explosion prevention
- The recognized rules applicable when connecting sterilized device components to non-sterile equipment!

⚠ WARNING

Danger of injury from electric shock!

To prevent any risk of electric shock, this device may only be connected to a 100–240 VAC power supply system having a protective conductor.

NOTICE

The SonicWeld Rx[®] system requires special precautions regarding electromagnetic compatibility (EMC), so when setting it up and putting it into operation, be sure to follow the EMC instructions provided in this user manual, see section 9 "Guidelines and Manufacturer's Declaration Regarding Electromagnetic Compatibility (EMC)", page 71.

NOTICE

Portable and mobile RF communication equipment can adversely affect the SonicWeld Rx[®] system.

4.3.2 Setting up SonicWeld Rx[®]

NOTICE

The SonicWeld Rx[®] ultrasonic unit (1, front fold-out page):

- must always be placed on a level, dry and solid surface.
 - must be set up and operated in the non-sterile area of the operating environment.
-
- Connect the handpiece (8, front fold-out page) to the connecting socket (2, front fold-out page). After inserting the plug of the connecting cable fully into the socket and rotate to lock it in place.

- The connecting cable of the handpiece is approximately 2.95 m long. If this length is not sufficient, an additional handpiece with a 6-m connecting cable is available as well.
- Screw the sonotrode manually in place on the handpiece and use the open-end wrench (12, front fold-out page) to check it for secure attachment (torque: max. 0.3 Nm). Note that the wrench is shipped non-sterile, so it needs to be sterilized prior to use.

NOTICE**Disconnection from the supply system!**

The SonicWeld Rx® must be set up so that it can be easily and reliably disconnected from the power supply.

- The ultrasonic unit features 2 connecting sockets for handpieces (2, front fold-out page). This enables you to connect 2 handpieces in order to work alternately with 2 sonotrodes (e.g. 2 standard sonotrodes or 1 standard sonotrode and 1 smoothing sonotrode).
- Plug the mains cable into the ultrasonic unit (5, front fold-out page) and then into a mains socket-outlet with ground contact. As soon as the unit has been connected to the power supply, it is automatically set to standby mode. Therefore, full disconnection is possible only by pulling the plug of the mains cable out of the socket-outlet.

⚠ WARNING**Explosion hazard!**

- The SonicWeld Rx® is not suitable for use in atmospheres enriched with flammable mixtures of anesthetics and air, oxygen or nitrous oxide.

NOTICE

- Make sure that no objects or appliances will be deposited on the ultrasonic unit.
- Observe the recognized rules when connecting sterilized system components to non-sterile equipment!

⚠ WARNING**Danger of injury due to malfunction!**

SonicWeld Rx® must not be placed close to or stacked with other devices (always maintain separation distance!), see section 9 "Guidelines and Manufacturer's Declaration Regarding Electromagnetic Compatibility (EMC)", page 71.

If operation in close proximity to or in a stack with other devices is unavoidable, the SonicWeld Rx® must be monitored to verify its proper performance in such a setup.

4.3.3 Changing the Display Language

The display language has been factory-preset to English. Users can select a different language by proceeding as follows:

- Place the ultrasonic unit (1, front fold-out page) on a level, dry and solid surface.
- Connect a handpiece (8, front fold-out page) to a connecting socket (2, front fold-out page), inserting the plug of the connecting cable fully into the socket.
- Connect the mains cable first to the unit (5, front fold-out page) and then to a grounding socket-outlet. Once connected to a mains socket-outlet, the unit automatically goes into standby mode.
- With the unit in standby mode, simultaneously operate the hand activation pushbutton (10, front fold-out page) and the On/Standby switch (4, front fold-out page) and hold them down until the display lights up.
- Using the **hand activation** pushbutton, you can now select among the following languages: German, English, French, Spanish, and Portuguese.
- Confirm your preferred language by operating the On/Standby switch. This turns off the unit automatically.
- Upon turning on the unit again with the On/Standby switch, the display language will have been changed as desired for future use.

4.3.4 Visual Check

- Check visually that all components used are free of defects.
 - Damaged components must not be used. Replace!

NOTICE

The tip of the sonotrode must not show any signs of mechanical damage. It must not be bent or rammed (typical damage after the handpiece has fallen down).

4.3.5 Automatic Functional Test


CAUTION

Danger of burns!

If the sonotrode is allowed to make skin contact during the functional test, this can cause burns.

Avoid contact with skin, eyes, etc.!

Self-test: Upon turning on the unit with the On/Standby switch (4, front fold-out page), the handpieces (if connected) are ready for a self-test. This is indicated for each handpiece by the flashing activation LED (9, front fold-out page) and the following symbol flashing on the corresponding side of the display:

The self-test for this handpiece will be started upon operating the hand activation pushbutton (10, front fold-out page) for the first time. Hold down the button until the display message changes from **Testing** to **Cooling**, the activation LED stops flashing and the acoustic signal changes. Be sure to keep the tip of the sonotrode out of contact with objects during this process. If the test is successful, the unit is automatically set to working mode.

- Checkmark ✓ displayed (checkmark right or left, depending on which handpiece has been tested): Means “no malfunction detected”.
- Generator-related error messages are displayed in text form, while handpiece states are indicated by pictogram on the respective side (depending on handpiece), see section 6.1 “Troubleshooting”, page 64.

4.4 Intraoperative Use

4.4.1 Sterile Working

NOTICE

Non-sterile delivery of the SonicWeld Rx[®] system

All components of the SonicWeld Rx[®] system are delivered non-sterile.

Sonotrodes, open-end wrench and handpieces with connecting cables are located in the sterile area of the operating environment, which is why they must be used in sterile condition.

- Make sure that the handpieces are dry.
- Ensure that the sonotrodes, the open-end wrench and the handpieces with connecting cables have been properly sterilized; see section 5 "Processing, Cleaning, Disinfection and Sterilization", page 60.

▲WARNING

Danger of infection by non-sterile handling!

Improper sterilization and non-sterile handling of the SonicWeld Rx[®] system can pose serious health risks to patients.

Make sure that the SonicWeld Rx[®] system components to be used are sterile, see section 5 "Processing, Cleaning, Disinfection and Sterilization", page 60.

The handpieces with connection cables, the sonotrodes and the open-end wrench are steam-sterilizable, see section 5 "Processing, Cleaning, Disinfection and Sterilization", page 60.

4.4.2 Preparatory Measures

NOTICE**Verify the system's readiness for operation prior to each use!**

Be sure to check prior to each intraoperative use that the SonicWeld Rx® is ready for operation!

CAUTION**Danger of burns!**

If the activation button (10, front fold-out page) has been operated, the sonotrode can cause burns when allowed to make tissue contact.

- Avoid contact with skin, eyes, etc.!
-

NOTICE**Glove damage by sonotrode tip**

The sonotrode must be handled with special care to avoid glove perforation and user injury by the tip of the sonotrode.

- Once the mains cable has been connected first to the unit and then to a mains socket-outlet with ground contact, the unit automatically goes into standby mode.
- Turn on the unit with the On/Standby switch (4, front fold-out page).
- Press the handpiece to start the functional test. Take care to keep the tip of the sonotrode out of touch with objects during this process!
- Upon determining the depth and diameter of the pilot hole to be drilled in accordance with the indication, select the appropriate SonicPin and drill bit. It must be possible to insert the bit up to the STOP mark.

The geometry of the drill bit and the design of the SonicPin have been perfectly matched to each other. The working end of the bit is designed so that it allows you to drill through the Resorb x® or Resorb xG mesh, plate, foil or membrane implants. The bit geometries also reflect the diameter of these implants.

We recommend selecting a drill bit that is at least 1 mm longer than the pin to be inserted.

- Hold the handpiece like a pencil, making sure that the activation pushbutton (10, front fold-out page) can be easily operated with the index finger.

4.4.3 Application

▲WARNING

Danger of infection in case of sterile glove damage!

There is a risk of sterile glove perforation when touching the sonotrode.

- When working with the sonotrode, be sure to avoid:
 - perforation of the sterile gloves
 - user or patient injury by contact with the hot tip of the sonotrode
 - To pick SonicPins from the clip magazine, insert the tip of the sonotrode into the SonicPin and withdraw the pin.
 - As the SonicPins attach magnetically to the tip of the sonotrode, they can easily be taken out of the clip magazine.
-


- Place the SonicPin attaching to the sonotrode tip into the borehole of the plate or mesh.
- Exert slight pressure on the SonicPin axially to the drilled hole.

- Operate the activation pushbutton (10, front fold-out page) until the SonicPin has been inserted into the hole to the intended depth under slight pressure, then stop. The blue indicator (9, front fold-out page) lights up on the handpiece during activation and an alert is sounded in addition.

Avoid any sonotrode tip movement for at least 3 s after inserting the SonicPin.

- Once the SonicPin has solidified, remove the sonotrode tip with a twist.


Resorb x® or Resorb xG mesh or foils can also be fixed as follows:

- First, implant a sufficient number of SonicPins in the bone as described above, using the standard sonotrodes (11, 13, front fold-out page) for pin insertion.

- Put the Resorb x® or Resorb xG membrane or foil in place over the inserted pins.
- Use the smoothing sonotrode (14, 15, front fold-out page) or the angled working end of the standard sonotrode (13, front fold-out page) to weld the Resorb x® or Resorb xG membrane or foil in place on the inserted pins.
 - Place the tip of the smoothing sonotrode onto the membrane/foil at the point of contact between the membrane/foil and the SonicPin underneath.
 - Exert just a little pressure on the membrane/foil.
 - Operate the activation pushbutton until the membrane liquefies visibly. As long as the pushbutton is pressed down, the blue activation LED (9, front fold-out page) lights up on the handpiece and an alert is sounded in addition.
 - Do not move the tip of the sonotrode for at least 3 s after operating the pushbutton.
 - Once the membrane/foil has solidified again, withdraw the tip of the sonotrode.
- Repeat this procedure to fix the Resorb x® or Resorb xG membrane or foil to all the remaining SonicPins.

The smoothing sonotrodes can also be used for smoothing out sharp or projecting edges of the membrane or foil.

- Place the tip of the smoothing sonotrode at the edge of the membrane or foil.
- Exert just a little pressure on the edge of the membrane or foil.
- Operate the activation pushbutton until the membrane liquefies visibly. As long as the pushbutton is pressed down, the activation LED (9, front fold-out page) lights up on the handpiece and an acoustic signal is heard in addition.
- Run the tip of the sonotrode very slowly along the edge of the membrane/foil during activation of the system.

NOTICE

If activation is continued for 60 s, the handpiece will be switched off automatically and the activation LED (9, front fold-out page) will go out.

4.5 Putting the Unit out of Operation

To put the unit fully out of operation:

- Operate the On/Standby switch (4, front fold-out page).
 - This sets the unit to standby mode and the display (3, front fold-out page) turns dark.
- Disconnect all system components (sonotrode from handpiece, handpiece with connecting cable from ultrasonic unit and mains cable from ultrasonic unit).

For all-pole disconnection of the unit from the power supply, pull the mains plug.

5 Processing, Cleaning, Disinfection and Sterilization

The warranty regarding handpiece cleaning and sterilization is limited to 250 reprocessing and sterilization cycles. Once the maximum number of cycles has been reached, a corresponding message will be displayed.

The recommendations of the Working Group Instrument Reprocessing (AKI) provide further information on how to properly reprocess reusable medical devices:

www.a-k-i.org > **AKI brochures** > **Red brochure**.

Further information on the reprocessing of medical devices can be found at:

www.klsmartin.com > **Service** > **Maintenance of Instruments**.

NOTICE

The responsibility for proper cleaning, disinfection and sterilization of the handpiece rests with the operator/user. Be sure to observe your national/local regulations, including potential restrictions.

Use appropriate and approved sterile packaging (e.g. in accordance with EN 868, ISO 11607) for sterilization, subsequent transportation and storage.

5.1 Cleaning/Disinfecting

5.1.1 Handpiece, by Machine

The handpiece of the SonicWeld Rx® is suitable for machine processing/thermal disinfection. It can be processed with the same programs that have been released for surgical instruments and implants. As regards cleaning, be sure to follow the instructions provided by the manufacturer of your washer disinfectant(s) (W/D) as well as those provided by the manufacturers of the cleaners and disinfectants used. The process (including loading) must guarantee sufficient removal of residues.

Only mildly alkaline cleaning agents are authorized for use.

- Remove the sonotrode from the handpiece.
- Disconnect the connecting cable with handpiece from the ultrasonic unit.
- Coil the connecting cable with handpiece in a circle having a diameter of at least 20 cm (17, rear fold-out page).
- We recommend using our tray basket for storing the connecting cable in the machine, see section 7 "Accessories & Spare Parts", page 69.
- Arrange the handpiece so inside the washer disinfectant that the sonotrode side of the handpiece points downward. This prevents any accumulation of rinsing liquid in the handpiece.

5.1.2 Handpiece, by Hand

- Remove the sonotrode from the handpiece.
- Disconnect the connecting cable with handpiece from the ultrasonic unit.

Clean the handpiece with a clean, lint-free cloth moistened with a commercially available disinfectant approved for use with instruments and based on ethanol (50/50) or isopropanol (70/30) (e.g. neodisher[®] MediClean by Dr. Weigert).

NOTICE

Risk of damage due to improper handling!

The handpiece must not be immersed in a disinfecting or ultrasonic bath.

Do not use any acetone-containing disinfectants for cleaning and disinfecting the handpiece.

5.1.3 Device, by Hand

Clean the SonicWeld Rx[®] ultrasonic unit with a clean, lint-free cloth moistened with a commercially available disinfectant approved for surface disinfection and based on ethanol or methanol (e.g. neoform MED AF by Dr. Weigert).

⚠ WARNING

Danger of death by electric shock!

- Ingress of liquids into the device must be avoided under any circumstances!
 - Before cleaning or disinfecting the SonicWeld Rx[®], the device must be fully disconnected from the power supply in all poles.
 - Should moisture have penetrated into the device, dry it after disconnecting it from the supply.
-

5.2 Visual Check for Damage

- Perform a visual check of all components before using the system again and before sterilizing the handpiece.
In particular, check the cable connections and verify that the unit and the sonotrode are undamaged.

▲WARNING**Danger of infection and injury by damaged components!**

- Damaged components must not be used!
-

5.3 Sterilization

▲WARNING**Danger of infection due to non-sterile handling!**

Improper sterilization and non-sterile handling of the SonicWeld Rx® system can lead to serious health hazards for patients.

Sterilization must be carried out according to a validated steam sterilization process, for example in a sterilizer complying with EN 285:2009 and ANSI/AAMI ST79, validated in accordance with ISO 17665-1:2006.

The fractionated (pulsing) vacuum method requires sterilization at 134°C (273.2°F)/2 bar with a minimum holding time of 5 min.

ANSI/AAMI ST79 recommends a minimum cycle time of 4 min at 132°C (270°F) for dynamic air removal steam sterilization cycles. Please follow the instructions of the user manual of your steam sterilizer.

- **Open-end wrench**

The wrench (12, front fold-out page) can be reprocessed in the same way as surgical instruments. A wrench in sterile packaging is required for each use.

- **Sonotrode**

The sonotrodes must be reprocessed in the same way as surgical instruments. Each reuse requires the sonotrodes to be available in new sterile packaging.

- **Handpieces with connecting cables**

The handpiece with a 2.95-m connecting cable (7, 8, front fold-out page) and the handpiece with a 6-m connecting cable are steam-sterilizable.

We recommend using our tray basket or, for small parts, the tray basket with cover, see section 7 "Accessories & Spare Parts", page 69.

NOTICE

Danger of patient infection or injury!

The warranty regarding cleaning and sterilization of the handpiece with connecting cable is limited to 250 reprocessing and sterilization cycles. Once the maximum number of cycles has been reached, a corresponding message will be displayed on the ultrasonic unit.

5.4 Maintenance

The maintenance measures required for the SonicWeld Rx® system are described in section 11 "Safety Checks (SC)", page 75.


6 Maintenance

6.1 Troubleshooting

6.1.1 Displays

Message displayed	Cause	Remedy
Fehler Stromversorgung [Power supply error] Temperaturmessung defekt [Temperature meter defective] Übertemperatur Generator [Generator overtemperature] Gerätefehler [Internal fault]	Generator error or fault.	Disconnect the unit from the power supply, and then reconnect it. If the problem persists, return unit to manufacturer, see section 1.1, page 45.
Störung Ultraschall [Ultrasound trouble]	Generator problem (during welding process).	Acknowledge by pressing the On/Standby switch. If the problem persists, disconnect the unit from the supply.
1. Handstück entlasten! [Relieve pressure on handpiece!] 2. On/Standby drücken! [Press On/Standby!]	System overload (sonotrode tip stuck/jammed, excessive pressure).	Acknowledge by pressing the On/Standby switch, then clear the handpiece problem and operate the activation pushbutton again
Schweißzeit > 60s [Welding time]	Maximum welding time of 60 s exceeded.	After a cooling-down time of 5 s, the unit is ready for use again.
1. Sonotrode festziehen! [Tighten sonotrode!] 2. On/Standby drücken! [Press On/Standby!]	Sonotrode seated too loosely.	Tighten sonotrode with wrench (12, front fold-out page), and acknowledge error.
Systemfehler [Internal fault]	Monitoring of all important program modules.	Disconnect the unit from the power supply, and then reconnect it. If the problem persists, return unit to manufacturer, see section 1.1, page 45.

6.1.2 Handpiece Statuses

Message displayed	Cause	Remedy
	No handpiece connected.	Connect handpiece.
	Temperature values measured on handpiece not plausible.	Replace handpiece.
	Max. temperature reached on handpiece.	Message disappears automatically as soon as the temperature drops below 70°C (158°F). Alternative: Replace handpiece.
	The handpiece cycle counter pictogram appears during self-test. (240 ≤ cycle counter < 250)	Message disappears automatically once the self-test has been completed.
	The flashing handpiece cycle counter pictogram appears during self-test. (cycle counter ≥ 250)	Following completion of the self-test, the message must be acknowledged with the On/Standby switch.
	Handpiece pushbutton defective.	Replace handpiece.

6.1.3 Potential Malfunctions

Trouble	Cause	Identification	Remedy
“Screeching” noises heard during the self-test or during use of the SonicWeld Rx® system.	Residual moisture or condensate on the contact points.	Clearly audible.	This does not compromise the proper functioning of the SonicWeld Rx® system. Noise will subside automatically during normal use of the device. It may be helpful to use longer drying times after sterilization.
Sonotrode tip will not detach from SonicPin. SonicPin fails to attach to sonotrode tip.	Bent or rammed sonotrode tip.	Visual check.	Replace sonotrode tip.
Sonotrode tip sticks to SonicPin.	Ultrasound applied too long.	Pin is generally deformed, sonotrode cannot be removed residue-free.	Observe required cooling-down time. Rotate the sonotrode axially (twist) to break the connection, see section 4.4, page 55.
	Bent sonotrode tip.	Visual check.	Replace sonotrode tip.
Handpiece indicator does not light up during activation.	LED life exceeded.	Visual check.	Replace the handpiece. A defective indicator on the handpiece has no adverse effect on ultrasound application.
Display dark.	No power supply.	Visual check.	Connect the system to the power supply and operate the On/Standby switch.
	Mains fuse defective.	Electrical test.	Replace mains fuse.
	Unit not switched on.	Visual check.	Press On/Standby switch until display lights up.

6.2 Replacements

6.2.1 Replacing the Sonotrode

- Use the wrench (12, front fold-out page) to remove the sonotrode (11, front fold-out page) from the handpiece by counterclockwise rotation; subsequently tighten the new sonotrode with the wrench (clockwise rotation).

6.2.2 Replacing the Mains Fuse

- Set SonicWeld Rx® to standby mode by pressing the On/Standby switch.
 - The unit can now no longer be activated and the display turns dark upon entering standby mode.
- Remove mains cable from mains socket-outlet.
- Pull out the fuse holder (6, front fold-out page).
- Remove the defective mains fuse from the fuse holder.
- Insert new mains fuse into the fuse holder, see section 7 "Accessories & Spare Parts", page 69, into the fuse holder.
- Push fuse holder back into the unit, making sure that it locks in place.
- Check the unit for proper functioning.

6.2.3 Replacing the Activation LED

The activation LED (9, front fold-out page) cannot be replaced by the user. However, a defective LED has no adverse effect on the proper functioning of the SonicWeld Rx® system as a whole.

6.2.4 Replacing the Mains Cable

- Set SonicWeld Rx® to standby mode by pressing the On/Standby switch. The unit can now no longer be activated.
 - The display turns dark upon entering standby mode.
- Unplug the mains cable first from the mains socket-outlet and then from the IEC inlet connector (5, front fold-out page).
- Use new country-specific mains cable (see section 7 "Accessories & Spare Parts", page 69) and connect it first to the unit (5, front fold-out page) and then to a grounding outlet.
- Check the unit for proper functioning.

6.2.5 Replacing the Handpiece with Connecting Cable

- Set SonicWeld Rx® to standby mode by pressing the On/Standby switch. The unit can now no longer be activated.
 - The display turns dark upon entering standby mode.
 - Remove the plug of the handpiece connecting cable from the handpiece connecting socket (2, front fold-out page).
 - Connect a new handpiece with connecting cable (7, 8, front fold-out page; see section 7 "Accessories & Spare Parts", page 69) to the SonicWeld Rx® system: insert the plug of the connecting cable fully into the handpiece connecting socket (2, front fold-out page) and rotate to lock it in place.
 - Operate the On/Standby switch to set the SonicWeld Rx® to working mode.
 - Check the unit for proper functioning.
-

 WARNING**Danger of injury in the event of using non-approved accessories!**

Be sure to use only the genuine parts available from the manufacturer as specified in section 7 "Accessories & Spare Parts", page 69. Using third-party accessories or spare parts can increase interference emission and reduce the interference immunity of the SonicWeld Rx®.

7 Accessories & Spare Parts

Designation	REF
Ultrasonic Unit SonicWeld Rx®	52-500-21-04
Handpiece with connecting cable (2.95 m)	52-500-23-04
Handpiece with connecting cable (6 m)	52-500-24-04
Standard sonotrode, straight	52-501-21-04
Standard sonotrode, angled	52-501-22-04
Smoothing sonotrode, straight	52-501-23-04
Smoothing sonotrode, angled	52-501-24-04
Open-end wrench for sonotrodes	52-502-01-04
Instructions for Use, SonicWeld Rx®	90-267-58-30
Mains fuse (T 2 AH, 250 V)	52-503-24-04
Mains cable Germany, France, Austria, Scandinavia	08-024-00-39
Mains cable Great Britain, Ireland, UAE	08-024-00-40
Mains cable Italy	08-024-00-30
Mains cable Japan	08-024-00-31
Mains cable Canada, Korea, Taiwan, USA	08-024-00-29
Mains cable Switzerland	08-024-00-28
Mains cable China	08-024-00-41
Storage facilities	
Tray basket for MiniSet Container System 277 x 171 x 54 mm (L x W x H) with cover	55-015-30-01
Tray basket	55-969-44-04
Tray basket cover	55-963-51-04
Storage and sterilization container	
MicroStop® – MiniSet Container System	55-861-70-04
Coding label with text	55-864-01-04
Logistics frame, red, for container	55-864-12-04
2 coding labels with text, for MiniSet Container System	55-864-05-04

8 Technical Data

Designation	SonicWeld Rx®	
Mains fuse (2 pieces)	T 2 AH, 250 V	
Power supply frequency	50–60 Hz	
Power supply voltage	100–240 V alternating voltage with protective ground contact	
Power input	70 VA	
Maximum power	17 W	
Operating frequency	30 kHz	
Protection class	I acc. to IEC 60601-1	
Degree of protection	IP42	
Suitability for use in atmospheres enriched with oxygen	The SonicWeld Rx® is not suitable for use in atmospheres enriched with flammable mixtures of anesthetics and air, oxygen or nitrous oxide.	
Operating mode	Continuous duty with intermittent loading 2 s operation 15 s break 5 repetitions 4 min. break	
Applied part	Handpiece with type BF sonotrode	
Dimensions	Width	260 mm
	Height	255 mm
	Depth	46 mm
Length of handpiece cable	2.95 m	
Length of long handpiece cable	6 m	
Number of reprocessing cycles for the handpiece	250	
Environmental conditions for transport and storage	Ambient temperature	-20°C to +45°C (-4°F to +113°F)
	Relative humidity	10–100%
	Atmospheric pressure	500–1,060 hPa
Environmental conditions for operation	Ambient temperature	+10°C to +40°C (+50°F to +104°F)
	Relative humidity	30–75%
	Atmospheric pressure	700–1,060 hPa
	Altitude	max. 4,000 m
	Conforms to 93/42/EEC	

9 Guidelines and Manufacturer’s Declaration Regarding Electromagnetic Compatibility (EMC)

Guidance and manufacturer’s declaration – Electromagnetic interference emission		
<p>The SonicWeld Rx® is intended for use in an electromagnetic environment as specified below. The user should assure that it is used in such an environment.</p>		
Emission test	Compliance level	Electromagnetic environment – Guidance
RF emissions according to EN 55011	Group 2	The SonicWeld Rx® uses electromagnetic energy exclusively for the performance of its internal functions. This may cause interference in nearby electronic equipment.
RF emissions according to EN 55011	Class B	The SonicWeld Rx® is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions according to IEC 61000-3-2	Complies	
Voltage fluctuations/flicker emissions according to IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – Electromagnetic interference emission

The SonicWeld Rx® is intended for use in an electromagnetic environment as specified below. The user should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guidance
Electrostatic discharge (ESD) according to IEC 61000-4-2	±6 kV contact discharge ±8 kV air discharge	±6 kV contact discharge ±8 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Fast electrical transients/bursts according to IEC 61000-4-4	±2 kV for mains power lines ±1 kV for input and output lines	±2 kV for mains power lines ±1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge according to IEC 61000-4-5	±1 kV differential-mode voltage ±2 kV common-mode voltage	±1 kV differential-mode voltage ±2 kV common-mode voltage	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short power interruptions and supply voltage fluctuations according to IEC 61000-4-11	<5% U_T (>95% dip in U_T) for ½ period 40% U_T (60% dip in U_T) for 5 periods 70% U_T (30% dip in U_T) for 25 periods <5% U_T (>95% dip in U_T) for 5 seconds	5% U_T (95% dip in U_T) for ½ period 40% U_T (60% dip in U_T) for 5 periods 70% U_T (30% dip in U_T) for 25 periods 5% U_T (95% dip in U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user requires continued operation during mains power interruptions, it is recommended to power the device from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field according to IEC 61000-4-8	3 A/m	3 A/m	Power-frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note	U_T is the a. c. mains voltage prior to application of the test level.		

Guidance and manufacturer's declaration – Electromagnetic interference emission

The SonicWeld Rx® is intended for use in an electromagnetic environment as specified below. The user should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – Guidance
Conducted RF according to IEC 61000-4-6 Radiated RF according to IEC 61000-4-3	3 V _{eff} 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V _{eff} 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	Portable and mobile RF communications equipment should be used no closer to any part of the device (including cables) than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.17\sqrt{P}$ $d = 1.17\sqrt{P}$ for 80 MHz to 800 MHz $d = 2.33\sqrt{P}$ for 800 MHz to 2.5 GHz where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol. 
Note 1	At 80 MHz and 800 MHz, the higher frequency range applies.		
Note 2	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio stations and AM and FM radio and TV stations, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level specified above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device. ^b Over the frequency range of 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF telecommunications equipment and the Gebrüder Martin SonicWeld Rx®

The SonicWeld Rx® is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF telecommunications equipment (transmitters) and the device as recommended below, depending on the maximum output power of the communications equipment.

Rated maximum output <i>P</i> of transmitter in watts	Separation distance according to frequency of transmitter (in m)		
	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.


Note 1	At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
Note 2	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

10 Disposal of the Device

Be sure to observe all national/local regulations and guidelines when disposing of the packaging material and potentially infectious items (e.g. removed metal components).

In designing the device, we tried to avoid using composite materials wherever possible. This allows a high degree of recycling. We therefore offer to take the device back for proper disposal and recycling.

Be sure to observe your national/local rules and regulations governing disposal!

	<p>Marking of electric and electronic equipment in accordance with Directive 2002/96/EC (WEEE Directive) and the German Electrical and Electronic Equipment Act (ElektroG)</p> <p>This symbol on the product or its packaging indicates that the product may not be disposed of as normal household garbage.</p>
---	---

11 Safety Checks (SC)

NOTICE

All safety checks must be documented in the medical device logbook together with the inspection and test results!

Once a year a safety check must be carried out acc. to IEC 62353 in order to verify that the unit duly conforms to safety regulations. The results of each check must be documented and the proper working condition of the unit confirmed by the person performing the check.

For frequency and scope of the safety check, see IEC 62353.

The points listed here are additional check points.

Frequency:	annually
Device class (MDD):	IIb
Inv. number:	

Device:	SonicWeld Rx®
Serial number:	
Operating hours:	

Visual inspection

OK **not OK**

- External surfaces of the unit (e.g. display, housing, handpiece)
- Mains cable
- Handpiece connecting cable
- Plugs
- Equipment complete
- Instructions for Use available
- All screws of the housing in place

<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>

Functional test

OK **not OK**

- On/Standby switch
- Display
- Activation pushbutton
- Activation LED during activation

<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>

Electrical safety (DIN EN 60601)

OK **not OK**

- Protective conductor resistance: $\leq 0.2 \Omega$
- Leakage current: $\leq 0.5 \text{ mA}$

<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>

Measures taken

- Technical Service contacted
- Device withdrawn from service

Date: _____

Date: _____

*Person performing
the checks/tests:*

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Rudolf Buck GmbH
KLS Martin UK Ltd.
KLS Martin Australia Pty Limited