Gebrauchsanweisung
Operating Instructions
Instrucciones de uso
Mode d'emploi
Istruzioni d'uso

Operationsleuchten	ML 1001/702/701/502/501
Operating Lights	ML 1001/702/701/502/501
Lámparas de quirófano	ML 1001/702/701/502/501
Luminaire opératoires	ML 1001/702/701/502/501
Lampade scialitiche	ML 1001/702/701/502/501

© by **Martin Medizin-Technik**

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© by **Martin Medizin-Technik**

This document is under the copyright protection of **Martin Medizin-Technik**.

This document must not be modified, extended or reproduced or transmitted to a third party, without prior written permission by **Martin Medizin-Technik**.

Martin Medizin-Technik

Ludwigstaler Str. 132
D-78532 Tuttlingen

Telephone: + 49-74 61/70 60
Telefax: + 49-74 61/70 61 93
E-mail: info@martin-med.com
Internet: www.martin-med.com

In case of technical questions or replacement parts are required, please contact our hotline service.

Martin Medizin-Technik

Hotline:
Telephone: + 49-74 61/70 63 43
Telefax: + 49-74 61/70 61 90
E-mail: msc@martin-med.com

Replacement parts:
Telephone: + 49-74 61/70 61 99
Telefax: + 49-74 61/70 61 90

Hotline service available on workdays from 8:00 a.m. to 5:00 p.m.

We reserve the right to make changes in designs and devices.

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1 Introduction

1.1 Foreword

In their construction and technical overall conception in the field of lighting engineering, the **Martin** Operating Lights* outstandingly fulfill the high requirements which are set by today's operating technology for the lighting of the operating area.

The lights bear the EC-sign, i.e. they meet the basic requirements stipulated in the EC-directive on medical products.

In addition, the lights ML 1001, ML 701 and ML 501 are UL approved.


These instructions for use will provide detailed information about how to put the operating lights in service and use them for the intended purpose.

Being a user or a service technician authorized by the manufacturer, please read these instructions carefully prior to using the light in the operating room, in order to make yourself acquainted with the handling and the operation of the operating light.

The manufacturer warrants that this operating lights have been properly and carefully tested before leaving the factory premises. In addition, we are prepared to assume continued responsibility for the safety, reliability and proper functioning of the operating lights as long as all checks and repairs required are carried out by technical staff that has been duly authorized by Martin to perform such work, and in as much as the operating lights are used only for the intended purposes.

Martin will assume no liability whatsoever for any damage that has been caused by making improper use of the product!

* The terms "operating light" or "light" respectively used in these operating instructions not only refer to the luminous element, but also to the mounting frame, the tracking and spring arms, the stands, the control boxes, etc.!

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1.2 Symbols within these
Operating Instructions

Texts marked by this symbol and set in bold letters or italics contain very important information, like for example, the prevention of hazards to health!

Please observe these texts under any condition!

Texts marked by this symbol and set in bold letters or italics contain very important information, like for example, the prevention of material damage!

Please observe these texts under any condition!

Texts marked by this symbol and set in bold letters or italics contain comments/ hints and other useful information.

1.3 General Information

1.3.1 Validity

The descriptions within these operating instructions refer to the **Martin** operating lights with the designation

- ML 1001 E/R;
- ML 702 HX *vario*;
- ML 701 E/R;
- ML 502 X *vario*;
- ML 501 E/R and

their respective versions and options.

If a description applies to all models of a given lamp type, no specific designation (e.g. "HX *vario*") will be indicated in the documentation.

Note the designation on the light, which is imprinted on the light head, and look for its specific description on the respective pages of these instructions.

The other chapters are valid for all lights.

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1.3.2 Proper Use for the Intended Purposes

The lights may only be used for purposes being in accordance with their design:

i.e., for lighting an operating or examination area by surgeons, doctors, or appropriately trained paramedics.

1.3.3 Entry Inspection

Please check the light(s) and the accessories for shipping damage and completeness immediately after having received them.

1.3.4 Complaints

Damage may only be claimed for if the vendor or forwarding company are immediately being informed.

Immediately establish a certificate of damage and send it to the nearest agent of the manufacturer or directly to the manufacturer.

If possible, use the original packages when returning goods.

Enclose the following accompanying documents with the returned goods:

- name and address of the owner
- type and device number
- description of the defect

1.3.5 Warranty

Our Standard Terms and Conditions effective at the time shall apply.

Within the warranty period, our competent service agencies or the manufacturer remedy all defects which were demonstrably caused by manufacturing faults or faulty material for free.

Parts subject to wear like lamps and replaceable hand grips are excluded from this guarantee.

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2 Safety Instructions

2.1 Safety Instructions for the Lights
ML 1001/702/701/502/501

All Martin operating lights are high-quality products that have been designed and manufactured in accordance with technical standards. They fully meet, at the time when they leave the factory, all safety-related requirements.

In order to preserve this condition, you as the user are requested to refrain or avoid all actions that may have a negative influence on the safety of the lights, and to have repaired damaged parts immediately or notify the manufacturer of the damage!

Damaged lights may be a life- or health-threatening hazard for the user and/or the patient!

For this reason never use damaged lights!

The lights are equipped with filter systems which prevent the radiation heat from the lights entering the operation area. Without these filters the lights may not be operated. Defective filter glass has to be replaced immediately.

When set to the maximum brightness and the smallest diameter the fields of light from light combinations should not be superimposed. There is a risk that too much heat would build up in the light field.

Only handle, switch on or off the lights in accordance with the information provided in these instructions for use!

The mounting and the repair of the lights may be carried out either by the manufacturer themselves or by a person or company explicitly authorized by the manufacturer!

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

In case the repair is executed by a person or company authorized by the manufacturer, the operator of the light is ask to demand a certificate containing the type and the extent of the repair by the person or company carrying out the repair.

This must contain

- date of the service,
- type and extent of the service,
- name of the company, and
- the signature of the service technician.

If repair work has been carried out by a third party, the equipment or parts repaired must be labeled with the repairer's identification mark.

After having carried out the mounting and service tasks, the service technician has to check all functions of the lights, spring and tracking arms, and the power supply units for their perfect functioning and operational safety!

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Laser pilot (Option)

The built-in laser pilot is a Class 2 laser that meets all requirements set by EN 60825/IEC 825/VDE 837. These requirements as well as pertinent regulations for the prevention of accidents (VBG 93, as far as Germany is concerned) must be observed in using the laser pilot.

The illustrated warning labels are located on the control panel where the laser pilot is activated.

The laser beam exit opening is located on the central handle. Whenever the operating light is repositioned with the help of this handle (or, in the case of the ML 1001, by using the hand rail), the laser beam will be switched on automatically. After repositioning has been completed, the laser pilot light will persist for another 2 s.

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Notes concerning prevention of accidents:

- *Due to the eyelid closure reflex, the laser radiation poses no danger to non-narcotized persons.*
- *In the case of narcotized patients, however, take care that the laser radiation will not hit the opened eye for more than 0.25 s. Cover the patient's eyes by using suitable means.*

The user has to assume full responsibility for the observation of the regulations mentioned above.

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2.2 Supplementary Safety Instructions for the Lights ML 702 and ML 502

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In contrast to the other lights, the ML 702 and ML 502 lights are not operated with low voltage, but with 300 V DC! This voltage is generated in the light head.

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In order to actuate the lamp in the light, up to 5000 V AC are generated in the light head and conducted via the focus of the lamp!

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In order to avoid a life-threatening electrical shock

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- *never dismantle, repair or modify any parts of these lights!*
- *never reach into the light head when the focus has been removed!*

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In case the light or any of its parts have been damaged,

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- *cut off the power supply of the light!*
- *secure the light or the parts against persons who might touch the light, if necessary, by putting up appropriate danger signs!*
- *immediately call the service technician!*

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For technical reasons, the ML 502 lights do not incorporate a reserve lamp. To avoid a complete loss of light,

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- *these lights should only be operated in combination with other lights!*
- *at least two lights should be directed at the surgical field!*

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3 Versions and Combinations of the Lights

3.1 Versions of the Lights

The lights are being delivered in different versions. These versions differ

- in terms of the **diameter of the light head**, e.g.:
 - ML 1001 corresponds to a light head diameter of 1000 mm
 - ML 701 corresponds to a light head diameter of 700 mm,
 etc.
- in terms of the **handling**, e.g.:
 - The version ML 701 E contains an electronical control system with a control panel at the hand grip of the light.
 - The version ML 701 R is a manual version with a control panel at the light head. As an option version R can be equipped with a separate brightness controller.
- in terms of the **type of mounting**, e.g.:
 - ceiling mounting
 - wall mounting
- in terms of the **options** by means of which the lights may be equipped, like for example:
 - laser pilot
 - medTV *pro plus*-system
 - external brightness controller
- in terms of the **lamp technology**, e.g.:
 - halogen lamp
 - metal halide lamp

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- The operating lights with the designation

- ML 1001
- ML 701
- ML 501

are equipped with one or two halogen lamp(s) respectively.

- The operating lights with the designation

- ML 702

are equipped with a halogen lamp and a metal halide lamp.

- The operating lights with the designation

- ML 502

are equipped with a metal halide lamp.

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The lights with the designation

- ML 1001
- ML 701
- ML 501

may be mounted as individual light but also as double, triple or fourfold combination (with identical or different equipment); the combination has no impact on the handling of the lights.

Lights of the type

- ML 702

can be used both individually or in combination with other lights (i.e. with additional ML 702 or ML 502 lights).

Lights of the type

- ML 502

should only be used in combination with other lights (see above).

The handling of the individual light is light-specific.

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4 Description of the Lights

By using the new **Martin** filter technology, these lights produce the **Martin** Super Cold Light.

4.1 Description of the Light ML 1001

The **Martin** ML 1001 operating light is available in two versions:

- The basic model ML 1001 E, where all functions can be electronically controlled via a central control panel in the tube connection piece.
- The ML 1001 R is the manual version with a control panel in the tube connection piece.

The **Martin** ML 1001 operating light provides the following accessory equipment:

Martin Anti-Drift System (ADS)

In the case of the **Martin** ML 1001, the operating light is arrested and kept in position by means of the built-in magnetic brakes. As soon as the user touches either the hand rail or the handle of the operating light, the magnetic brakes located in the rotary joints are released and the operating light can be repositioned with ease and precision.

Upon releasing the rail or the central handle, the operating light will be locked securely in its new position by the **Martin** Anti-Drift System.

Martin Laser pilot (only ML 1001 E)

The laser pilot built into the operating light exactly marks the center of the light field. This makes it particularly easy for the user to focus the light field with utmost precision even if the operating area is very small. Whenever the operating light is moved/repositioned, the laser pilot will be activated. When the operating light has come to a standstill in the course of repositioning, the laser pilot continues to be in operation for another two seconds.

The laser pilot can be deactivated via the control panel.

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Martin medTV^{pro plus}/medTFT^{pro}

The **Martin medTV^{pro plus}** is an imaging system intended for documentation and training purposes.

The high-performance camera may be mounted either on the operating light itself or on a separate camera support.

For image display, the optional **Martin medTV^{pro}** flat-panel monitor system is available.

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4.2 Description of the Light ML 702

The **Martin ML 702** operating light only comes in one available type. It may either be used individually or in combination with other ML 702 or ML 502 lights.

By applying a new kind of lamp technology, i.e. a metal halide lamp, a high light efficiency is achieved while there is only a minor heat build-up in the wound area. Compared to the ordinary operating lights with halogen lamp technology, the costs for operation and service are thus relatively low.

The characteristic technical features of this lamp are:

- A light efficiency of about 45 % of the absorbed energy.
- Light intensity of up to 160,000 lux with a power input of 70 watt.
- A very low emission of infrared rays.
- With a service life of more than 5,000 hours, a five to ten times longer service life than that of halogen lamps.

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Start-up behavior

Upon switching on the operating light, the halogen lamp is activated and moved into central position by the electronic system. At the same time, the metal halide lamp is ignited. In the case of a cold start, the metal halide lamp requires approx. 2 minutes to achieve 90% of its light output. After this start-up phase, the halogen lamp is switched off automatically and the metal halide lamp is moved into central position (which takes approx. 3 s).

In case of power failure, less than 0.5 seconds (according to the current standard) the change-over board automatically switches to battery power supply, but the electric arc is fully maintained (no intensity fluctuations).

If the light is switched off and immediately switched on again, the metal halide lamp requires a warm-start period of approx. 10 minutes before it can be re-ignited by the ballast. During this warm-start phase, the halogen lamp will provide the light output (just like during a cold start). The yellow indicator light on the control unit will indicate this.

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**4.3 Description of the Light
ML 701**

The **Martin** ML 701 operating light is available in two versions:

- The model ML 701 E, where all functions can be electronically controlled via a central control panel located on the handle of the operating light.
- The ML 701 R is the manual version with a control panel located at the light head. As an option when used as a ceiling light version R can be equipped with a separate brightness controller.

Optionally, the ML 701 can be combined with the following accessory equipment:

Martin Laser pilot (only ML 701 E)

The laser pilot built into the operating light exactly marks the center of the light field. This makes it particularly easy for the user to focus the light field with utmost precision even if the operating area is very small. Whenever the operating light is moved/repositioned, the laser pilot will be activated. When the operating light has come to a standstill in the course of repositioning, the laser pilot continues to be in operation for another two seconds. The laser pilot can be deactivated via the control panel.

Martin medTV pro plus/medTFT pro

The **Martin** medTV pro plus is an imaging system intended for documentation and training purposes. The high-performance camera may be mounted either on the operating light itself or on a separate camera support. For image display, the optional **Martin** medTV pro flat-panel monitor system is available.

**External brightness control
(only ML 701 R ceiling light)**

The brightness of the lights can be varied with the external brightness control. The controller is, e.g. installed in a wall operating panel.

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Start-up behavior

As soon as the light is switched on, the metal halide lamp will be ignited automatically. In the case of a cold start, this bulb requires approximately two minutes to reach 90 % of its light output and be ready for use.

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If the light is switched off and immediately switched on again, the metal halide lamp requires a warm-start period of approx. 10 minutes before it can be re-ignited by the ballast.

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**4.5 Description of the Light
ML 501**

The **Martin** ML 501 operating light is available in two versions:

- Model ML 501 R is the manual version with an operating panel located at the light head. As an option when used as a ceiling light version R can be equipped with a separate brightness controller.
- The ML 501 E, where all functions can be electronically controlled via a central control panel located at the light head.

Optionally, the ML 501 can be combined with the following accessory equipment:

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Martin medTV *pro plus*/medTFT *pro*

The **Martin** medTV *pro plus* is an imaging system intended for documentation and training purposes. The high-performance camera may be mounted either on the operating light itself or on a separate camera support. For image display, the optional **Martin** medTV *pro* flat-panel monitor system is available.

**External brightness control
(only ML 501 R ceiling light)**

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4.6 Options for the Lights ML 1001/702/701/502/501

These operating lights have optionally, at the ceiling tube, a card for the commutation mains/battery.

A double LED on the soffet indicates optically the power source actually in use.

Green LED: Power supply (mains)

Red LED: Emergency power supply (battery)

Note:

If neither of the LED's is lit, the battery does not supply any power.

The LEDs are connected to the emergency relais with a three-pin-plug. For light combinations, the LED does always concern the main light.

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5 Operating the Lights

5.1 Operating the Lights
ML 1001 E

All functions of the **Martin** ML 1001 E operating light can be easily controlled via a separate control panel that has been arranged in a straightforward manner. The various functions available have been marked with a symbol. To activate a function via the respective pushbutton, no more than a touch of the finger is required.

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- 1 Operating light OFF/ON
(or use wall switch, as the case may be)
- 2 Laser pilot (option)
Laser stand-by button
(the green LED will come on)
- 3 Malfunction indicator lamps:
 - Main lamp I defective
 - Main lamp II defective
- 4 Decreasing light intensity
- 5 Light intensity indicator
- 6 Increasing light intensity
- 7 Reducing light field diameter
- 8 Light field indicator
- 9 Increasing light field diameter

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
The modern **Martin** operating lights are equipped with microprocessor-controlled electronics. In case an error occurs, one of the LEDs in the luminous intensity indicator (5) will be flashing. At the same time, either an LED display moves back and forth within the light field indicator (8) or all LEDs are flashing at the same time.

In case of under- or overvoltage, the two malfunction indications (3) will be flashing.

In all the cases mentioned above please contact the technical service since the repair work required may only be carried out by service technicians.

Note:

In general, this does not lead to a failure of the lights. It is not necessary to switch off the lights and the operation may be terminated. In certain cases, the nominal luminous intensity will be reduced to 70 % approx.

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
Martin Anti-Drift System (ADS)

The **Martin** Anti-Drift System (ADS) allows easy and high-precision positioning of the operating light. As soon as the operating light has come to a standstill, it will be fixed in position by the built-in magnetic brakes. These brakes will be automatically released whenever the hand rail surrounding the operating light (or the central handle) is touched. Thus, repositioning is possible with utmost ease.

Accordingly, the operating light will again be arrested in its new position as soon as the user withdraws his hand from the rail or handle.

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Martin Laser pilot

To put the laser pilot into operation, press the pushbutton 2 on the control panel. As soon as the central handle or the hand rail is touched, the laser will be activated.

Upon releasing the rail or handle, the laser light will persist for another 2 seconds.

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**5.2 Operating the Lights
ML 1001 R**

All functions of the **Martin** ML 1001 R Operating Light can be easily controlled via a control panel located in the tube connection piece. The various functions available have been marked with a symbol. To activate a function via the respective pushbutton, no more than a touch of the finger is required.

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- 1** Operating light OFF
(or use wall switch, as the case
may be)
- 2** Operating light ON
(first switch on the wall switch
if provided)
- 3** Indicator lamp for monitoring the
two light sources.
If one of the main light sources has
become defective, this indicator
lamp will light up (red).

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Note:

The **Martin** ML 1001 R also features the **Martin Anti-Drift System (ADS)**. For a detailed description of this system, please refer to Section 5.1.

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**5.3 Operating the Lights
ML 702 and ML 502**

All functions of the **Martin** ML 702/502 Operating Light can be easily controlled via a control panel located in the tube connection piece. The various functions available have been marked with a symbol. To activate a function via the respective pushbutton, no more than a touch of the finger is required.

The symbols on the panel represent the following functions:

- 1 Operating light OFF
(or use wall switch, as the case
may be)
- 2 Operating light ON
(first switch on the wall switch
if provided)
- 3 Main light source monitoring
indicator (only ML 702).
This indicator lights up (yellow)
whenever the halogen lamp is
active. In normal operation, this is
the case during the start-up phase
of the metal halide lamp (cold start:
approx. 2 min; warm start: approx.
10 min). However, the indicator also
lights up in the event of a metal
halide lamp failure.
- 4 Indicator signifying normal
operation. It will show up green
once the metal halide lamp has
been ignited and the light is ready
for use.
- 5 Increasing light intensity
- 6 Indicator (yellow): lights up when
the "increase light intensity" button
has been pressed.
- 7 Indicator (yellow): lights up when
the "decrease light intensity" button
has been pressed.
- 8 Decreasing light intensity

Please note:

If the indicator (3) light remains on,
it signals that both the metal halide
lamp and the halogen lamp should
be replaced as soon as possible.

To increase or decrease the light
field diameter, rotate the sterilizable
handle (1) clockwise or anticlockwise
as appropriate.

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**5.4 Operating the Lights
ML 701 E and ML 501 E**

All functions of the **Martin** ML 701 E/ML 501 E operating light can be easily controlled via a separate control panel that has been arranged in a straightforward manner. The various functions available have been marked with a symbol. To activate a function via the respective pushbutton, no more than a touch of the finger is required.

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- 1 Operating light OFF/ON
(or use wall switch, as the case
may be)
- 2 Laser pilot (option)
Laser stand-by button
(the green LED will come on)
- 3 Malfunction indicator lamps:
- Main lamp I defective
- Main lamp II defective
- 4 Decreasing light intensity
- 5 Light intensity indicator
- 6 Increasing light intensity
- 7 Reducing light field diameter
- 8 Light field indicator
- 9 Increasing light field diameter

Please note:

If one of the malfunction indicator lamps (3) has come alight, the respective halogen lamp should be replaced as soon as possible.

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5.5 Operating the Lights ML 701 R and ML 501 R

All functions of the **Martin** ML 701 R/ML 501 R Operating Lights can be easily controlled via a separate control panel that has been arranged in a straightforward manner. The various functions available have been marked with a symbol. To activate a function via the respective pushbutton, no more than a touch of the finger is required.

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- 1 Operating light OFF
(or use wall switch, as the case
may be)
- 2 Operating light ON
(first switch on the wall switch
if provided)
- 3 Indicator lamp for monitoring the
two light sources.
If one of the main light sources has
become defective, this indicator
lamp will light up (red).

Please note:

If the malfunction indicator lamp (3) has come alight, the respective halogen lamp should be replaced as soon as possible.

While the halogen lamps have a rated life of between 600 and 1000 burning hours, this average service life may be significantly reduced in consequence of mechanical vibrations.

To increase or decrease the light field diameter, rotate the sterilizable handle (1) clockwise or anticlockwise as appropriate.

When used as a ceiling light as an option the **Martin** ML 701 R/ML 501 R operation light can be equipped with an external brightness control which can be used to vary the brightness of the light between 35 % and 100 %. In accordance with the design an additional switch which is used to switch the light on and off is located next to the controller.

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5.6 Operating the Wall-Mounted Operating Lights ML 501 E and ML 501 R


All wall-mounted operating lights are provided with a wall-mounting control box containing a power switch. This switch must be set to ON prior to operating the operating light via the control panel located at the light head of the operating light.

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The wall-mounting control box also contains the line fuses.

	220 - 240 V	100 - 125 V
ML 501	2,0 AT	4,0 AT

	220 - 240 V	100 - 125 V
ML 501	2.0 AT	4.0 AT

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5.7 Replacing the Sterilizable Handle

To remove the sterilizable handle, just press the locking bar located on the upper end of the handle to release the lock. The handle can then be with-drawn with ease.

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When inserting the handle, be sure that it clicks in place audibly. To achieve this, it might be necessary to rotate the handle just a little clockwise or anticlockwise.

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ML 1001	
Order No.	Description
86-910-02	Halogen lamp, 24 V/250 W (6 pcs.)
86-900-00	Sterilizable handle for operating lights not equip- ped with a laser pilot
86-990-04	Sterilizable handle for operating lights equip- ped with a laser pilot

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ML 702	
Order No.	Description
86-907-02	Lamp replacement kit (1 pcs.)
86-900-00	Sterilizable handle for operating lights not equip- ped with a laser pilot

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ML 701	
Order No.	Description
86-907-01	Halogen lamp, 24 V/150 W (6 pcs.)
86-900-00	Sterilizable handle for operating lights not equipped with a laser pilot
86-990-04	Sterilizable handle for operating lights equipped with a laser pilot

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7 **Maintenance and Servicing of the Light**

To keep the operating light in good working order, all service and maintenance work must be carried out with due care.

7.1 **Cleaning and disinfection of the externally accessible parts**

All parts of the lights which can be reached from the outside can be cleaned using a neutral cleaning agent based on tensides (handwashing liquid, neutral cleaning agent).

To disinfect the lights neutral disinfection agents (e.g. surface disinfectants) based on aldehydes, quaternary ammonium compounds, etc. can be used.

In order to prevent humidity from penetrating into the light and the tracking and spring arms, do **not** spray these parts but only wipe them with cloths moistened with detergents or disinfectants.

7.1.1 **Cleaning, disinfection and sterilisation of the handles**

See Ch. 5.7 for information regarding the assembly / disassembly of the handle for sterilisation.

The replaceable handle can be treated mechanically with the same programmes approved for surgical instruments. The replaceable handle is thermally stable and can be thermally disinfected at 93 °C.

The manufacturer's instructions for the cleaning equipment, as well as those of the manufacturer of the treatment agent have to be observed.

In order to prevent spots or the buildup of deposits we recommend that fully demineralised water is used for the final rinse.

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Cleaning can be carried out with neutral or alkaline cleaning agents. Cleaning agents which contain active chlorine should be avoided due to the fact that reactions between the chlorine and plastic can be expected.

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The handle can be sterilised with steam. The sterilisation is carried out in a fractional vacuum process at 134 °C / 5 minutes or 121 °C / 20 minutes (e.g. in a steriliser which complies with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993 and has been validated according to EN 554/ISO 13683).

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During the sterilisation the handle may not be mechanically burdened as this could lead to permanent deformation.

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Warning:
Under no circumstances may the handle be treated in a hot air steriliser. This would result in damage!

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Service life:
The replaceable handles are subject to natural wear and tear if sterilized frequently. If fatigue occurs (e.g., cracks forming; discoloration), the handle must be replaced.

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7.2 Maintenance

At regular intervals all functions of the operating light need to be checked for proper functioning and operating reliability.

Every 24 months a Safety Check **must** be carried out by an authorized person. For details, please refer to Section 11.

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7.2.1 Replacing the Halogen Lamps of the Lights ML 1001/701/501

The procedure for replacing a defective halogen lamp is identical for the ML 1001, ML 701 and ML 501 operating lights.

WARNING!

Never perform work on the operating light while it is switched on!

First, manually remove the three countersunk screws on the focusing unit, then withdraw the focusing unit completely from the light housing.

The electronic lamp monitoring system ensures that the defective lamp is automatically moved aside with the help of a servo-motor provided for this purpose. At the same time, the second light source is moved into the focus of the reflector.

Accordingly, always replace the defective lamp located on the edge of the lamp holder.

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Before removing the (multi-blade) lamp diaphragm (1), please check the position of the lamp's connecting lines.

Unlock the diaphragm (1) rotating it counterclockwise through 45 degrees, then withdraw it in upward direction.

Remove the halogen lamp (2) and the metal halide lamp (3), disposing of them in accordance with national and regional regulations.

Note:
Never touch the glass bodies of the new halogen lamps with your bare hands!
When holding a new lamp, always grasp it on the connector or base!

Remove the lamp's protective cover, then plug the halogen lamp (using the lamp insertion aid) and the metal halide lamp into the appropriate lamp socket.

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
Warning:

Whenever a halogen lamp is replaced, the lamp sockets must be checked for potential defects. Remember that lamp sockets are wearing parts!

When replacing a lamp, always check the cables for scorches in the lamp socket area. If such defects are found, please call your serviceman.

Install the diaphragm by placing it onto the support and rotating it clockwise by approx. 45 degrees. Be sure that the connector (6) is parallel to the surface (4) of the diaphragm attachment.

Take care that the lamp's connecting cables are installed properly. These lines should be run in loops and kept as close to the guide bar (5) as possible, irrespective of the lamp carrier's position (retracted or extended). Be sure to prevent the connecting cables from rubbing against the filter glass unit when the lamp carrier moves.

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When mounting the focusing unit, be sure to insert the plug (1) into the rectangular cutout (2) of the assembly support.

To do this, manually insert the focusing unit in such a manner that it is flush with the assembly support. Then lock the unit in place by rotating the three countersunk-head screws clockwise.

Please switch on the operating light to check it for proper functioning. When cold-starting the ML 702, the indicator light (see Section 5.3, pos. 3 on the control panel) must go off after approx. 2 minutes. At the same time, the metal halide lamp is moved into central position.

For further information on the operating light's start-up behavior, please refer to Section 4.2.

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**7.2.3 Replacing the Metal Halide
Lamp of the Light ML 502**

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Remove the metal halide lamp (1), by pulling out the lamp base (2), then dispose of it in accordance with applicable regional/national regulations.

Note:

Only use genuine Martin metal halide lamps (see section 6, "Accessories").

Never touch the glass body of the new lamp with your bare hands!

When holding a new lamp, always grasp it on the connector or base!

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8 Trouble Shooting

8.1 Trouble Shooting ML 1001/701/501

Symptom	Potential Cause	Remedy
Operating light cannot be switched on.	Wall-mounted switch not set to ON.	Set wall-mounted switch to ON.
Malfunction indicator lamps I and II alight on the control panel.	Both halogen lamps defective.	Replace both halogen lamps.
Malfunction indicator lamp I or II alight on the control panel.	One halogen lamp defective.	Replace the halogen lamp.
Laser pilot beam not visible.	Laser pilot stand-by button (on control panel) not activated. Sterilizable handle without glass insert on the handle.	Press laser pilot stand-by button Use sterilizable handle with glass insert.
Brightness control not working.	Focusing unit not properly inserted.	Insert focusing unit properly.
Light field diameter cannot be adjusted.	Focusing unit not properly inserted.	Insert focusing unit properly.

8.2 Trouble Shooting ML 702/ML 502

Symptom	Potential Cause	Remedy
- Light of the operating light is too dark or too yellow or blue	- The lamp in the light has not yet attained the operating temperature	- Wait until the operating temperature has been reached (cold start: approx. 2 min; warm start: approx. 10 min)
	- End of the lamp service life attained	- Mind yellow indicator lamp (only ML 702) - Replace the metal halide lamp
	- Focus error	- Set the brightness/light field with the focus handle
- Yellow indicator lamp lights up (only ML 702)	- The lamp in the light has not yet attained the operating temperature	- Wait until the operating temperature has been reached (cold start: approx. 2 min; warm start: approx. 10 min)
	- Metal halide lamp defective	- Replace the metal halide lamp

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9 Specifications

	ML 1001	ML 702	ML 701	ML 502	ML 501
Voltage	110/230 V AC 50 - 60 Hz / 24 V DC				
Nominal output of the transformer	410 VA	300 VA	220 VA	168 VA	168 VA
Light source (lamps)	24 V/250 W	90 V/70 W*	24 V/150 W	90 V/70 W	24 V/120 W
Max. central illumination intensity E_C	145 klx	160 klx	140 klx	120 klx	90 klx
Total radiation intensity	514 W/m ²	550 W/m ²	497 W/m ²	412 W/m ²	315 W/m ²
Diameter of light field d_{10}	200 mm	190 mm	160 mm	180 mm	140 mm
Diameter of light field d_{50}	105 mm	95 mm	83 mm	90 mm	80 mm
Diameter of light field d_{10} small / large	200 / 350	190 / 270	160 / 250	180 / 240	140 / 200
Residual illumination intensity with a "Schatter"*** in reference to E_C	64 % (92800 lx)	40 % (64000 lx)	40 % (56000 lx)	10 % (12000 lx)	10 % (9000 lx)
Residual illumination intensity with two "Schattern"*** in reference to E_C	53 % (76850 lx)	41 % (65600 lx)	41 % (57400 lx)	41 % (49200 lx)	41 % (36900 lx)
Residual illumination intensity in a standard tube in reference to E_C	87 % (126150 lx)	100 % (160000 lx)	100 % (140000 lx)	100 % (120000 lx)	100 % (90000 lx)
Residual illumination intensity in a standard tube with a "Schatter"*** in reference to E_C	50 % (72500 lx)	40 % (64000 lx)	40 % (56000 lx)	10 % (12000 lx)	10 % (9000 lx)
Residual illumination intensity in a standard tube with two "Schattern"*** in reference to E_C	42 % (60900 lx)	41 % (65600 lx)	41 % (57400 lx)	41 % (49200 lx)	41 % (36900 lx)
Illumination depth L1 + L2	1240 mm	1420 mm	1180 mm	1300 mm	1200 mm
Colour temperature	4300 K				
Colour reproduction index R_a	93				
Lamp lifetime	600 - 1000 h***	approx. 5000 h***	600 - 1000 h***	approx. 5000 h***	600 - 1000 h***
Protection class 	Type B according to IEC 601-1				
Environmental operating conditions - Storage temperature - Operating temperature	- 10 °C through + 50 °C for an air humidity of 20 % through 90 % (not condensing) - 5 °C through + 40 °C for an air humidity of 30 % through 80 % (not condensing)				

* Metal halide lamp: 90 V / 70 W; halogen lamp: 24 V / 120 W.

** Standard-form (according to DIN EN 60601-2-41) which corresponds to the shadow produced by the surgeon's head.

*** Depending on mechanical vibration intensity, voltage fluctuations, and duty cycles.

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10 Meaning of the Serial Number

The serial number provided on the rating plate indicates the following:

1. Product
e.g.
ML: Martin operating light
2. Type
e.g.
1001 E: Light head 1,000 mm, with electronic light function control
3. Hardware version
e.g.
05: Hardware version 05
4. Software version
e.g.
00: Software version 00
5. Year of manufacture
e.g.
01: Built in 2001
6. Mounting type
e.g.
D: Ceiling-mounting
7. Consecutive numbering (serial number)
e.g.
105: Operating light no. 105
Each year, the serial number starts anew with "001".

In the event of inquiries or complaints, please always mention the complete serial number of the light.

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11 Regular Safety Checks

Overview of checks and tests to be carried out

At least every 24 months, the unit must be subjected to the following checks and tests (which need to be performed by fully qualified technicians who, owing to their training, expertise and practical experience, are in a position to carry out such work reliably and independently):

- Visually check the operating light and its accessories for mechanical damage that might impair the proper functioning of the unit. Also check that all safety labels are fully readable.
- Perform functional test in accordance with operating instructions.
- Carry out electrical checks in accordance with the test report sheet for periodic safety checks (see Service Manual).
- Measure the supply voltage at maximum brightness.

All safety checks and related results must be documented in an equipment log.

Martin operating lights may be checked only by Martin's own qualified service personnel or by specially authorised persons/firms!

Important note:
If the operating light is not functioning properly or not safe to operate, it must be repaired immediately by an authorized person! If this is not possible, the light must be withdrawn from service and suitable measures must be taken to prevent its use! At the same time, the operator of the light must be informed on these measures.

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Edizione Italiana copyright © - Tutti i diritti riservati.
Copyright by Gebrüder Martin GmbH & Co. KG
Gebrüder Martin GmbH & Co. KG